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ACCUVAC Rescue

Modified characteristics following PCB replacement

From serial number SN 25821 (and for a replacement in the course of a repair) the ACCUVAC Rescue device is supplied with a new electronics component. Once the electronics component has been replaced, some aspects of ACCUVAC Rescue will have modified characteristics. **These are not faults**, but characteristics of the new electronics component. The differences in the characteristics of ACCUVAC Rescue are explained below.

No.	Action performed	Old electronics component WM no. 10680, Index G	New electronics component WM no. 10680, Index H
1	Connecting cable disconnected from charging port	- Device switches off - O/I LED stays on for another 4 min.	- Device switches off - O/I LED goes out
2	A new battery has been fitted but not yet initialized	In pumping mode: - 10 % LED flashes During charging: - 10 % LED flashes If device switched off: - 10 % LED does not flash	In pumping mode: - 10 % LED flashes During charging: - 10 % LED does not flash, as it is incorporated in the "charging in progress" light If device switched off: - 10 % LED flashes for another 4 minutes
3	Excessive/inadequate temperature at battery	- Charging interrupted - "Charging in progress" light goes out	- Charging interrupted - "Charging in progress" light still on
4	Battery initialized in accordance with instructions for use: - press keys -0.05/-0.1/ -0.15 bar - Wait until the battery has discharged and pump switches off	10 % LED goes out	10 % LED remains on until the battery has subsequently been fully charged once
5	The battery is charged following initialization	During charging: - "Charging in progress" light on - Red LED on as long as battery capacity < 10 %	During charging: - "Charging in progress" light on - Red LED is incorporated in the "charging in progress" light and does not come on separately
6	The pump is operated at a specified pressure of -0.8 bar and with an open intake connector	- LED at -0.8 bar flashes - All other LEDs indicating pressure are off	- LED at -0.8 bar flashes - All other LEDs indicating pressure are off. The LED at -0.05 bar flickers occasionally, as the pump has a greater output and is thus already generating a slight vacuum with the intake connector open.
7	Battery being charged	Electronics do not make any noise	Individual electronics components make slight hissing noises. This is not a malfunction.

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Geräte für Medizin GmbH+Co. KG