

LIFE-BASE III MCI

Portable system for oxygen therapy in MCI and disasters

It is a major challenge for medical staff to provide oxygen to patients in the event of Mass Casualty Incidents (MCI) and disasters. First they have to give immediate care to large numbers of seriously injured patients. Support for their work comes from handy, portable and flexible inhalation units such as LIFE-BASE III MCI. The 2-liter oxygen cylinder integrated in the portable system bridges the gap until a complete external oxygen system with large cylinders can be set up.



LIFE-BASE III MCI,
WM 8090



LIFE-BASE III MCI,
WM 8257



LIFE-BASE III MCI,
WM 8226

As soon as the oxygen logistics system is ready for use, any LIFE-BASE III MCI can be connected to the external oxygen supplies by means of different quick-coupling systems and connection hoses. Designed with emergency conditions in mind, LIFE-BASE III MCI can simultaneously supply oxygen to as many as five patients quickly and efficiently.

- Simultaneous delivery of oxygen to several patients with a single system
- Simple inhalation can be expanded to ventilation function with use of MEDUMAT Easy ventilator
- Inhalation options via MODUL Oxygen or oxygen administrator
- Crash-resistant LIFE-BASE III MCI (as per EN 1789) fits in every vehicle equipped with the wall mounting BASE-STATION III or BASE-STATION 3 NG.
- Can be used in combination with other Weinmann products (e.g., with O₂ Box Rescue)
- Integrated gurney and bedrail holder
- Quick response with direct flow selection
- External oxygen feed possible
- DIN coupling for supply from external devices (WM 8257)

LIFE-BASE III MCI Portable System	WM-No.	NATO Stock Number
LIFE-BASE III MCI with MODUL Oxygen and 3-way administrator (pre-set to 10 l/min)	8226	–
LIFE-BASE III MCI with 4-way administrator with dosage setting stops	8257	–
LIFE-BASE III MCI with MEDUMAT Easy, MODUL Oxygen and 3-way administrator (pre-set to 10 l/min)	8090	6530-72-001-3397

1



2



3



Accessories

- 1 **Wall mounting BASE-STATION 3 NG to install LIFE-BASE III in rescue vehicles – WM 8190**
- 2 **3-meter hose for external gas supply with Walther sealing nipple and O₂ DIN connector – WM 22303**
- 3 **PEEP valve – WM 3215**

International connections available upon request. (bayonet (AFNOR) / AGA/BS 5682).

Our complete offering of system solutions for Emergency, Transport and Disaster Medicine can be found at www.weinmann.de

Technical Data LIFE-BASE III Portable System



Device dimensions (W x H x D in mm):	465 x 295 x 145	Protective bag:	removable	
Weight (empty, in kg):	3.0	Fulfills standards:	DIN EN 1789	
WM No.	Number of patients who can be cared for simultaneously	Number of patients who can be treated with inhalation pre-set at 10 l/min	Number of patients who can be treated with variable inhalation of 0 to 15 l/min	Number of patients who can be treated with demand flow or ventilation (via MEDUMAT Easy)
WM 8257	4	-	4	-
WM 8226	4	3	1	-
WM 8090	5	3	1	1

Upon request, other inhalation flows can be pre-set (e.g., one outlet at 4 l/min, one at 6 l/min, and one at 8 l/min)

LIFE-BASE III



WM 8257

LIFE-BASE III MCI with dosage setting stops on 4-way administrator

Consisting of:

- 4-way administrator with dosage setting stops (WM 3899) with:
 - Outlets: 4 x UNF 9/16 male thread, 0-15 l/min, dosage setting stops, incl. 4 x angled connection nozzles UNF 9/16 (WM 31130)
 - Inlet / Outlet: 1 x G 3/8 male thread, 1 x T-shaped piece without non-return valve Walther (WM 8187), 1 x O₂ coupling as per DIN 13260 (WM 30026)
- 1 x 2-liter oxygen cylinder (WM 1822)
- 1 x OXYWAY Fix III (WM 30301)
- 4 x inhalation masks with connection hose (WM 1429)
- 1 x pressure hose, 3 meters, sealing nipple (Walther) and O₂-DIN extraction connector to central gas system outlet
- LIFE-BASE III portable system, pre-assembled (WM 8041)
- Protective bag (WM 8295)



WM 8226

LIFE-BASE III MCI with MODUL Oxygen and 3-way administrator

Consisting of:

- MODUL Oxygen (WM 22360)
- 3-way administrator pre-set to 10 l/min per outlet (WM 3834)
- 2-liter oxygen cylinder (WM 1822)
- OXYWAY Fix III (WM 30301)
- Pressure hose (WM 5737)
- Pressure hose (WM 22320)
- Adapter for MODULs (WM 22265)
- 4 inhalation masks with connection hose (WM 1429)
- LIFE-BASE III portable system, pre-assembled (WM 8041)
- Protective bag (WM 8295)



WM 8090

LIFE-BASE III MCI with MEDUMAT Easy, MODUL Oxygen and 3-way administrator

Consisting of:

- MEDUMAT Easy (WM 28000)
- MODUL Oxygen (WM 22360)
- 3-way administrator pre-set to 10 l/min per outlet (WM 3834)
- 2-liter oxygen cylinder (WM 1822)
- OXYWAY Fix III (WM 30301)
- Administrator with Walther quick-coupling (WM 8220)
- Pressure hose (WM 22869)
- Adapter for MODULs (WM 22265)
- 4 inhalation masks with connection hose (WM 1429)
- LIFE-BASE III portable system, pre-assembled (WM 8041)
- Protective bag (WM 8295)

■ **Germany** ■ Weinmann Geräte für Medizin GmbH+Co.KG · P.O. Box 540268 · D-22502 Hamburg · KronsaaIsweg 40 · D-22525 Hamburg · E: info@weinmann.de · www.weinmann.de · T: +49-(0)40-5 47 02-0 F: +49-(0)40-5 47 02-461 ■ Center for Production, Logistics, Service · Weinmann Geräte für Medizin GmbH+Co.KG · Siebenstücken 14 · D-24558 Henstedt-Ulzburg · T: +49-(0)4193-88 91-0 F: +49-(0)4193-88 91-450 ■ **Switzerland** ■ G. Weinmann AG – Neuenhof · T: +41-(0)56-416 41 11 · E: info@ch.weinmann.de · www.weinmann.de ■ **France** ■ Weinmann S.A.S. – Paris-Igny T: +33-(0)1 69 35 53 20 · E: info@fr.weinmann.de · www.weinmann.de ■ **Russia** ■ Weinmann SPb GmbH – Saint-Petersburg · T: +7-(812)6 33 30 82 · E: info@ru.weinmann.de · www.weinmann.de ■ **China** ■ Weinmann GmbH+Co.KG (Shanghai Rep. Office) · T: +86 21 32 26 26 98 · E: info@cn.weinmann.de · www.weinmann.de

3140.1-01-EN-0909-1 © Weinmann, Hamburg. Duplication of any kind only with the express permission of Weinmann. All rights to design and specification modifications reserved.