

MEDUMAT Transport – Charging Options

Product to be charged	Power supply voltage	
	12 volt	100-240 volt
 <p>MEDUMAT Transport single devices WM 28300 and WM 28400</p>	 <p>WM 28356</p>	 <p>WM 28305</p>
 <p>MEDUMAT Transport on LIFE-BASE light WM 28340 and WM 28350</p>	 <p>WM 28356</p>	 <p>WM 28305</p>
 <p>MEDUMAT Transport on LIFE-BASE 1 NG WM 9625 and WM 9620</p>	 <p>WM 8214 or WM 9640 or WM 9107</p>	 <p>WM 9655 or WM 9117 each with WM 28305</p>
 <p>MEDUMAT Transport on LIFE-BASE 4 NG WM 9600 and WM 9605</p>	 <p>WM 8195 or WM 9107</p>	 <p>WM 8237 or WM 9117 each with WM 28305</p>
 <p>Akkupack PLUS WM 28385</p>	<p>no direct charge via 12-V vehicle's electrical system</p>	 <p>WM 28305</p>

3167.1-09-EN-1109-1 © Weinmann, Hamburg. Duplication of any kind only with the express permission of Weinmann. All rights to design and specification modifications reserved.



Germany ■ Weinmann Geräte für Medizin GmbH+Co.KG · P.O. Box 540268 · D-22502 Hamburg · Kronsaalsweg 40 · D-22525 Hamburg · E: info@weinmann.de · www.weinmann.de · T: +49-(0)40-5 47 02-0
 F: +49-(0)40-5 47 02-461 ■ Center for Production, Logistics, Service · Weinmann Geräte für Medizin GmbH+Co.KG · Siebenstücken 14 · D-24558 Henstedt-Ulzburg · T: +49-(0)4193-88 91-0
 F: +49-(0)4193-88 91-450 ■ **Switzerland** ■ G. Weinmann AG – Neuenhof · T: +41-(0)56-416 41 11 · E: info@ch.weinmann.de · www.weinmann.de ■ **France** ■ Weinmann S.A.S. – Paris-Igny
 T: +33-(0)1 69 35 53 20 · E: info@fr.weinmann.de · www.weinmann.de ■ **Russia** ■ Weinmann SPb GmbH – Saint-Petersburg · T: +7-(812)6 33 30 82 · E: info@ru.weinmann.de · www.weinmann.de
 ■ **China** ■ Weinmann GmbH+Co.KG (Shanghai Rep. Office) · T: +86 21 32 26 26 98 · E: info@cn.weinmann.de · www.weinmann.de

emergency | homecare | diagnostics

