



Sleep at last.

SOMNObalance e – auto-CPAP device with OPP technology for balanced therapy

- OPP (Obstructive Pressure Peak) to distinguish between central and obstructive apneas
- CPAP and auto-CPAP mode with optional softPAP pressure relief
- Epoch-based, needs-oriented pressure adjustment
- Small, light, quiet and pressure-constant (option: SOMNOaqua humidifier can be integrated)
- Patient-oriented information menu in display – success of therapy is visible
- Extensive statistical data for therapy monitoring
- Mask test function

www.weinmann.de ■ T: +49-(0)40-5 47 02-102