

PRESS RELEASE

Transport ventilation redefined

Date: 29 August 2008

MEDUMAT Transport for Emergency and Critical Care Transport Ventilation

Bremen – An innovative ventilator that can be used in emergency and critical care transport ventilation has been launched by Weinmann of Hamburg. MEDUMAT Transport, a new development in the proven MEDUMAT ventilation system, offers a multitude of additional features that make device use safer and simpler and ensure top-quality ventilation during critical care transport.

The versatile device, perfect for all areas of primary and secondary care, from transport of emergency, acute care and intensive care patients, has a unique integrated emergency mode, an innovative operating concept, a very large color TFT display, intuitive user guidance, several individually selectable types of ventilation and expanded ventilation monitoring. Combined with the tried and tested LIFE-BASE carrying systems, MEDUMAT Transport offers a broad range of uses in rescue services and critical care transport, trauma room and intra-hospital transfers.

Fast, safe and simple use of the emergency mode by even less-experienced operators is ensured with pre-settings for infants, children and adults. The user interface is similar in design to other MEDUMAT models. In emergency mode the device provides doctors and rescue personnel with all the information available in conventional emergency ventilators. In operating mode for differentiated ventilation, all settings are available for the types SPV (pressure-controlled) and SVV (volume-controlled). MEDUMAT Transport thus offers all the potential uses of a state-of-the-art hospital device. A total of eight different ventilation types can be selected: PCV, BiLevel, CPAP, ASB, SIMV, IPPV, SVV, SPV. The NIV option for non-invasive ventilation can be activated for all types of ventilation.

More information:

http://www.weinmann.de/en/home_emergency/beatmung_und_monitoring/medumat_transport/

(1799 characters)

Background: Weinmann

Weinmann Geräte für Medizin GmbH + Co. KG concentrates on its product lines Homecare, Emergency and Diagnostics. In these areas the Hamburg-based family business offers diagnostic, therapeutic and also life-saving devices and system solutions of the highest quality. Today it is in the hands of the Griefahn and Feldhahn families and its headquarters have been located in Hamburg-Stellingen since 1977.

Today the company has about 500 employees. With more than 100 patent families and 90 trademarks, Weinmann is among the leading providers of medical devices for home care, diagnostics and emergency medicine in Germany. Sales in 2007 were 65.5 million Euros. Weinmann does business in more than 50 countries and operates branch offices in France, Switzerland, Thailand, China, Russia, Australia and New Zealand.

Press Contact:

Weinmann Geräte für Medizin GmbH+Co.KG ■ Kronsaalsweg 40 ■ 22525 Hamburg ■ Germany

Frau Juliane Papendorf ■ Press Officer

T: +49 (0)40-5 47 02-562 ■ F: +49 (0)40-5 47 02-469

E: j.papendorf@weinmann.de ■ www.weinmann.de