

PRESS RELEASE

easyOX

Date: March 6, 2008

New Oxygen Concentrator by Weinmann

Hamburg – Patients with COPD, fibrosis or emphysema have to inhale additional oxygen for at least 16 hours per day. In the new oxygen concentrator easyOX, the Hamburg-based medical technology firm Weinmann offers a new device of proven Weinmann quality that combines maximum patient safety with a good price/performance ratio.

The easyOX is equipped with a Thomas compressor and is available in two different models. The single flow meter model is for use by one patient at home and the double meter model is for two patients in a hospital or nursing home. Patients receive oxygen in a continuous flow of one to five liters per minute, which is regulated by a flow meter integrated in the housing. A humidifier can be attached to the device's oxygen port if desired.

"For the safety of the patient, the easyOx is equipped with optical and acoustic alarms that monitor oxygen concentration, overpressure or low pressure and excess temperature," said Weinmann Product Manager Agim Imeri.

The device weighs 27.5 kg and runs quietly at a low sound level of 45 dB(A). The easyOX can be operated with 230 volts at 50 or 60 Hz and can therefore be used around the world.

More information at: http://www.weinmann.de/en/home_homecare/sauerstoffmedizin/easyox/

(1283 characters)

Background: Weinmann

Weinmann Geräte für Medizin GmbH + Co. KG concentrates on its product lines Homecare, Emergency and Diagnostics. In these areas the Hamburg-based family business offers diagnostic, therapeutic and also life-saving devices and system solutions of the highest quality. Today it is in the hands of the Griefahn and Feldhahn families and its headquarters have been located in Hamburg-Stellingen since 1977.

Today the company has about 500 employees. With more than 100 patent families and 90 trademarks, Weinmann is among the leading providers of medical devices for home care, diagnostics and emergency medicine in Germany. Sales in 2006 were 63.2 million Euros. Weinmann's Production, Logistics and Service Center was established in Henstedt-Ulzburg close to Hamburg in 2006. Weinmann does business in more than 50 countries and operates branch offices in France, Switzerland, Thailand, China, Australia, New Zealand and Venezuela.

Press Contact:

Weinmann Geräte für Medizin GmbH+Co.KG ■ Kronsaalsweg 40 ■ 22525 Hamburg ■ Germany

Juliane Papendorf ■ Press Officer

T: +49 (0)40-5 47 02-562 ■ F: +49 (0)40-5 47 02-469

E: j.papendorf@weinmann.de ■ www.weinmann.de