

PRESS RELEASE

Transport Ventilation redefined

Date: 14 May 2008

MEDUMAT Transport for Emergency and Transport Ventilation

Fulda – An innovative ventilator for use in emergency and hospital transport ventilation is being presented by the Hamburg-based medical technology company Weinmann at this year's RETTmobil. MEDUMAT Transport, the newest development in the proven MEDUMAT ventilation system, offers several extra features that make its use safer and simpler and guarantee high-quality ventilation of acute care patients.

The unique integrated emergency mode, innovative operating concept, very large color TFT display, intuitive user guidance and numerous individually selectable ventilation types and enhanced ventilation monitoring make the device the perfect multi-talent for use in primary and secondary care, from emergency to intensive-care transport. Combined with the proven LIFE-BASE carrying systems, MEDUMAT Transport offers a broad range of uses in rescue services and acute care transport, trauma room and intra-hospital transfers.

Fast, safe and simple use of the emergency mode by even less-experienced operators is ensured with pre-settings for infants, children and adults. The user interface is similar in design to other MEDUMAT models. In emergency mode the device provides doctors and rescue personnel with all the information available in conventional emergency ventilators. In operating mode for differentiated ventilation, all settings are available for the types SPV (pressure-controlled) and SVV (volume-controlled). MEDUMAT Transport thus offers all the potential uses of a state-of-the-art hospital device. A total of eight different ventilation types can be selected: PCV, BiLevel, CPAP, ASB, SIMV, IPPV, SVV, SPV. The NIV option for non-invasive ventilation can be activated for all types of ventilation.

Additional information:

http://www.weinmann.de/en/emergency_home/ventilation_and_monitoring/medumat_transport/

(1886 characters)

Background: Weinmann

Weinmann Geräte für Medizin GmbH + Co. KG concentrates on its product lines Homecare, Emergency and Diagnostics. In these areas the Hamburg-based family business offers diagnostic, therapeutic and also life-saving devices and system solutions of the highest quality. Today it is in the hands of the Griefahn and Feldhahn families and its headquarters have been located in Hamburg-Stellingen since 1977.

Today the company has about 500 employees. With more than 100 patent families and 90 trademarks, Weinmann is among the leading providers of medical devices for home care, diagnostics and emergency medicine in Germany. Sales in 2006 were 63.2 million Euros. Weinmann's Production, Logistics and Service Center was established in Henstedt-Ulzburg close to Hamburg in 2006. Weinmann does business in more than 50 countries and operates branch offices in France, Switzerland, Thailand, China, Australia, New Zealand and Venezuela.

Press Contact:

Weinmann Geräte für Medizin GmbH+Co.KG ■ Kronsaalsweg 40 ■ 22525 Hamburg ■ Germany

Juliane Papendorf ■ Press Officer

T: +49 (0)40-5 47 02-562 ■ F: +49 (0)40-5 47 02-469

E: j.papendorf@weinmann.de ■ www.weinmann.de