

PRESS RELEASE

easyOX

Date: 7.11.2006

New, efficient oxygen concentrator by Weinmann

(Hamburg) – Oxygen long-term therapy may only be successful, if applied on a daily basis for a minimum of 16 hours. It is therefore of utmost important to equip the patient with the most easy, safe and comfortable designed technology. Offering exactly that is easyOX, the new and highly efficient oxygen concentrator by Weinmann, Hamburg. The device comes in two models, made for export out of Germany only: as single Flowmeter for the patient at home, or as double flowmeter for two patients mainly for deployment in Hospitals or Care Facilities.

Two Flowmeter-tubes are fixed to the device and supply a oxygen flow of between 1 and 5 l/min. Already three minutes after having been put on easyOX reaches maximum oxygen concentration, which at 1 to 4 l/min is at 95 (+1/-3) Vol.-% of O₂ – guaranteeing highest therapeutic efficiency. The level of concentration is permanently watched over by an integrated O₂-Sensor. If required, an air humidifier may be attached to the O₂-outlet.

For the safety of the patient easyOX is equipped with optical and acoustic Alarms, controlling the O₂-concentration, over- and under-pressuring, as well as in case of excessive temperature. In addition, easyOX weighs only 27,5 kg and at 45 dB(A) is very quite. Thanks to its voltage compatibility of between 230 V at 50 or 60 Hz easyOX may be used all around the world. Weinmann offers a 3-year guarantee period.

(1380 Zeichen)

Weinmann Geräte für Medizin GmbH+Co. KG, Hamburg, offers high-quality system solutions in diagnostic, therapeutic, as well as life saving medical technology. The company concentrates on its business units Emergency, Homecare and Diagnostics.

Press Contact:

Weinmann Geräte für Medizin GmbH+Co. KG
Kronsaalsweg 40, 22525 Hamburg
Contact Person Press: Tobias Drewling
T: +49 (0)40 - 5 47 02 - 0
F: +49 (0)40 - 5 47 02 - 469
E: info@weinmann.de
www.weinmann.de