

PRESS RELEASE

O₂ Box Rescue

Date: 7.11.2006

New, mobile Oxygen-System for up to 2000 litres of supply: O₂ BOX Rescue

(Hamburg) – In case of a MIC (Mass Incident of Casualties) rescuers are depending on medical equipment being immediately ready, reliable, robust and easy to use. Weinmann's new O₂ BOX Rescue is made to meet exactly these prerequisites.

Moreover, O₂ Box Rescue allows the individual supply of oxygen for up to 4 patients at the same time. Upon opening, O₂ Box Rescue's manifold offers a variety of outlets, which may be, for example, used for a central Gassupply, as well as the appliance of an airpack, or demand-valve with ventilation bag - thereby allowing for maximum in flexibility. In addition, O₂ Box Rescue is the only mobile box with a supply of up to 2000 litres of O₂.

Using Walther-latches with the pressure hose and pressure reducer enables a speedy and easy change of O₂-bottles. For those bottles not in use the G 3/8" outlet and the coupling type Walther allow separate usages. National adapters and individually compiled mountings are available upon request.

(962 Zeichen)

Weinmann Geräte für Medizin GmbH+Co. KG, Hamburg, offers high-quality system solutions in diagnostic, therapeutic, as well as life saving medical technology. The company concentrates on ist business units Emergency, Homecare and Diagnostics.

Press Contact:

Weinmann Geräte für Medizin GmbH+Co. KG
Kronsaalsweg 40, 22525 Hamburg
Contact Person Press: Tobias Drewling
T: +49 (0)40 - 5 47 02 - 0
F: +49 (0)40 - 5 47 02 - 469
E: info@weinmann.de
www.weinmann.de