

# LIFE-BASE light

*The clever transport solution for all MEDUMAT ventilators*

Lightweight, small and compatible with all external oxygen sources: LIFE-BASE light, optimized for use in emergency medical services and hospitals. When LIFE-BASE light is connected to an external oxygen source, you have a complete system ready for immediate use.

**The lightweight portable system**



## Convenient transport and reliable operation

- extremely lightweight
- extremely small size requires less space
- compatible with all external oxygen sources
- varied fastening options for pole and standard rail mounting
- fulfills requirements of EN 1789 (not applicable for pole mounting)

Can be combined with:

- Weinmann's modules
- MEDUMAT Easy
- MEDUMAT Easy CPR
- MEDUMAT Standard
- MEDUMAT Standard a
- MEDUMAT Transport

## Overview of LIFE-BASE light Combination Options

Available Variants	Supply Schedule	
Portable system <b>LIFE-BASE light</b> with <b>MEDUMAT Easy CPR</b> and accessories bag <b>WM 28230</b>	<ul style="list-style-type: none"> <li>- LIFE-BASE light portable system</li> <li>- MEDUMAT Easy CPR</li> <li>- Accessories bag for LIFE-BASE light</li> </ul>	<ul style="list-style-type: none"> <li>- 3-m pressure hose with plug for central gas system outlet</li> <li>- Shoulder strap (optional)</li> </ul>
Portable system <b>LIFE-BASE light</b> with <b>MODUL CPAP</b> , small starter set for CPAP therapy and accessories bag <b>WM 28240</b>	<ul style="list-style-type: none"> <li>- LIFE-BASE light portable system</li> <li>- MODUL CPAP</li> <li>- Accessories bag for LIFE-BASE light</li> <li>- 3-m pressure hose with plug for central gas system outlet-Shoulder strap (optional)</li> </ul>	<ul style="list-style-type: none"> <li>- Small starter set for CPAP therapy, consisting of:               <ul style="list-style-type: none"> <li>- CPAP-Boussignac for MODUL CPAP</li> <li>- CPAP/NIV mask Size M</li> </ul> </li> <li>- Headgear</li> </ul>
Portable system <b>LIFE-BASE light</b> with <b>MEDUMAT Easy CPR</b> , <b>MODUL CPAP</b> , small starter set for CPAP therapy and protective bag with shoulder strap <b>WM 28235</b>	<ul style="list-style-type: none"> <li>- LIFE-BASE light portable system</li> <li>- MEDUMAT Easy CPR</li> <li>- MODUL CPAP</li> <li>- Protective bag for LIFE-BASE light with shoulder strap</li> <li>- 3-m pressure hose with plug for central gas system</li> </ul>	<ul style="list-style-type: none"> <li>- Small starter set for CPAP therapy, consisting of:               <ul style="list-style-type: none"> <li>- CPAP-Boussignac for MODUL CPAP</li> <li>- CPAP/NIV mask Size M</li> </ul> </li> <li>- Headgear</li> </ul>
Portable system <b>LIFE-BASE light</b> with <b>MEDUMAT Standard a</b> and protective bag with shoulder strap <b>WM 28245</b>	<ul style="list-style-type: none"> <li>- LIFE-BASE light portable system</li> <li>- MEDUMAT Standard a</li> <li>- MODUL CPAP</li> <li>- Protective bag for LIFE-BASE light with shoulder strap</li> </ul>	<ul style="list-style-type: none"> <li>- 3-m pressure hose with plug for central gas system outlet</li> <li>- Holder for patient valve</li> </ul>
Portable system <b>LIFE-BASE light</b> with <b>MEDUMAT Transport</b> <b>WM 28350</b> <b>WM 28340</b>	<ul style="list-style-type: none"> <li>- LIFE-BASE light portable system</li> <li>- MEDUMAT Transport with CO<sub>2</sub>-measurement (WM 28350)</li> </ul>	<ul style="list-style-type: none"> <li>- MEDUMAT Transport without CO<sub>2</sub>-measurement (WM 28340)</li> </ul>